2569769

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new case, keep base data

Version 2.7en
2012-12-03

1 Administrative information	
To which NCA(s) is this report being sent?	
Type of report	
<input type="radio"/> Initial report <input type="radio"/> Follow-up report <input checked="" type="radio"/> Final report	
Date of this report	
2023-07-14	
Reference number assigned by the manufacturer	
Absolute Pro LL May 23, 2022	
FSCA reference number assigned by NCA	
Incidence reference number assigned by NCA	
Name of the co-ordinating NCA Competent Authority (if applicable)	
ANSM	

2 Information on submitter of the report	
Status of submitter	
<input type="radio"/> Manufacturer <input type="radio"/> Authorised Representative within EEA and Switzerland <input checked="" type="radio"/> Others: (identify the role)	
Local Abbott Regulatory Contact	

3 Manufacturer information		new
Name		
Abbott Vascular Santa Clara		
Contact Name		
Lindsey Bell		
Address		
3200 Lakeside Drive		
Postcode	City	
95054	Santa Clara, CA	
Phone	Fax	
+1 951-914-5259	+1 951-914-3995	
E-mail	Country	
qahotline@abbott.com	US - USA	

4 Authorised Representative Information

new

Name Abbott Vascular International BVBA	
Contact Name Ani Karapetyan	
Address Park Lane, Culliganlaan 2B	
Postcode 1831	City Diegem
Phone + 32 2 714 1411	Fax +32 2714 1412
E-mail AVEMEAREgulatoryCompliance@av.abbott.com	Country BE - Belgium

5 National contact point information

new

National contact point name Abbott Laboratories LLC (AR)	
Name of the contact person Evgeniia Sharapova	
Address 125171, Russia, Moscow, Leningradskoe sh., 16A, b.1	
Postcode 125171	City Moscow
Phone +7 (965) 249-27-92	Fax +7 (495) 258-42-81
E-mail evgeniia.sharapova@abbott.com	Country RU

6 Medical device information

new

Class <input type="radio"/> AIMD Active implants <input type="radio"/> MDD Class III <input checked="" type="radio"/> MDD Class IIb <input type="radio"/> MDD Class IIa <input type="radio"/> MDD Class I <input type="radio"/> IVD Annex II List A <input type="radio"/> IVD Annex II List B <input type="radio"/> IVD Devices for self-testing <input type="radio"/> IVD General	
Nomenclature system (preferable GMDN) GMDN	Nomenclature code 47932
Nomenclature text Multiple peripheral artery stent, bare-metal	
Commercial name/ brand name / make Absolute Pro LL Peripheral Self-Expanding Stent Sy	
Model number All model numbers	Catalogue number
Serial number(s) Not Applicable	Lot/batch number(s) All lots
Device Mfr Date	Expiry date

Notified Body (NB) ID-number 2797
Accessories / associated devices (if applicable) Not Applicable
Software version number (if applicable) Not Applicable

7 Description of the FSCA

Background information and reason for the FSCA

Abbott is initiating a Field Safety Action for all lots of Absolute Pro LL Peripheral Self Expanding Stent System (PSESS). Abbott has confirmed reports of mechanical locking, stent deployment and partial stent deployment failures, resulting from unintended excessive force used to deploy the stent. To reduce occurrences of deployment failures and associated outcomes, Abbott is notifying all users of the potential causes and risks.

Description and justification of the action (corrective / preventive)

The Absolute Pro LL PSESS includes a self-expanding stent that is pre-mounted on an over-the-wire Delivery System. The Absolute Pro LL PSESS is intended for the stenting of peripheral arteries as an adjunct to percutaneous transluminal angioplasty (PTA) and for the palliation of malignant strictures in the biliary tree. Deployment related issues are occurring at a rate of 0.27%. Investigation of these events has identified particular use conditions that can increase the likelihood of a stent deployment failure or potentially worsen the outcomes should a partial stent deployment occur. In the event the stent is partially deployed, the inability to fully release or capture the stent has led to surgery or additional intervention. Potential patient effects include dissection / tissue damage, foreign body in the patient, and occlusion. In one case following a procedure, a series of cascading events that included surgery resulted in a patient's death. Therefore Abbott determined that field safety notice was required.

Advice on actions to be taken by the distributor and the user

Customers are advised to:

- Read through the Safety Notice
- Share the information with other personnel associated with Absolute Pro LL procedures in the organization
- If these products have further been distributed/transferred, notify those customers
- Sign and return the provided Effectiveness Check Form
- Report any occurrence of product performance issues or patient adverse events to Abbott

Progress of FSCA , together with reconciliation data (Mandatory for a Final FSCA)

2023/07/14

All the effectiveness check forms were signed and returned by the concerned customers.

Time schedule for the implementation of the different actions

2023/04/14 The IFU has been updated. All EU product built on or after June 1, 2023, will contain the updated IFU.

Attached please find

- Field Safety Notice (FSN) in English
- FSN in national language
- Others (please specify)

FSN Status

- Draft FSN
- Final FSN

The medical device has been distributed to the following countries:

within the EEA and Switzerland

- | | | | | | | | |
|--|--|--|--|--|--|--|--|
| <input checked="" type="checkbox"/> AT | <input checked="" type="checkbox"/> BE | <input type="checkbox"/> BG | <input checked="" type="checkbox"/> CH | <input checked="" type="checkbox"/> CY | <input checked="" type="checkbox"/> CZ | <input checked="" type="checkbox"/> DE | <input checked="" type="checkbox"/> DK |
| <input type="checkbox"/> EE | <input checked="" type="checkbox"/> ES | <input checked="" type="checkbox"/> FI | <input checked="" type="checkbox"/> FR | <input checked="" type="checkbox"/> GB | <input type="checkbox"/> GR | <input type="checkbox"/> HU | <input checked="" type="checkbox"/> IE |
| <input type="checkbox"/> IS | <input checked="" type="checkbox"/> IT | <input type="checkbox"/> LI | <input type="checkbox"/> LT | <input type="checkbox"/> LU | <input type="checkbox"/> LV | <input type="checkbox"/> MT | <input checked="" type="checkbox"/> NL |
| <input type="checkbox"/> NO | <input checked="" type="checkbox"/> PL | <input checked="" type="checkbox"/> PT | <input checked="" type="checkbox"/> RO | <input checked="" type="checkbox"/> SE | <input checked="" type="checkbox"/> SI | <input checked="" type="checkbox"/> SK | <input type="checkbox"/> TR |

Candidate Countries

- HR

All EEA, candidate countries and Switzerland

Others:

See attached list of additional countries where the device has been distributed.

8 Comments

Commercial name: Absolute Pro LL Peripheral Self-Expanding Stent System

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorised representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.

Signature

I affirm that the information given above is correct to the best of my knowledge

print

check

send XML-data by E-Mail